

## Product datasheet for **RP300022**

### **MIR124-1 Rat qPCR Primer Pair (MI0000893)**

#### **Product data:**

<b>Product Type:</b>	qPCR Primer Pairs
<b>Locus ID:</b>	100314291
<b>Forward Sequence:</b>	CGTGTTACAGCGGACCTTG
<b>Reverse Sequence:</b>	GAACATGTCTGCGTATCTC
<b>ACCN:</b>	MIMAT0004728
<b>Synonyms:</b>	MI0000893; MIMAT0004728; MIR124-1 rno-mir-124-1; rno-miR-124-5p
<b>Components:</b>	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
<b>Quality Control:</b>	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
<b>Storage:</b>	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)