

## Product datasheet for **RP300071**

### MIR181b-1 Rat qPCR Primer Pair (MI0000926)

#### Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100314040
Forward Sequence:	AACATTCATTGCTGTCGGTGG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000859
Synonyms:	MI0000926; MIMAT0000859; MIR181b-1 rno-mir-181b-1; rno-miR-181b-5p
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)