

Product datasheet for **RP300072**

MIR181b-2 Rat qPCR Primer Pair (MI0000927)

Product data:

| | |
|-------------------|--|
| Product Type: | qPCR Primer Pairs |
| Locus ID: | 100314159 |
| Forward Sequence: | AACATTCATTGCTGTCGGTGG |
| Reverse Sequence: | GAACATGTCTGCGTATCTC |
| ACCN: | MIMAT0000859 |
| Synonyms: | MI0000927; MIMAT0000859; MIR181b-2 rno-mir-181b-2; rno-miR-181b-5p |
| Components: | 1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions) |
| Quality Control: | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage: | The primer mix is stable for one year from date of shipping. Store at -20°C. |



[View online »](#)