

Product datasheet for **RP300126**

MIR218a-2 Rat qPCR Primer Pair (MI0000957)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100314057
Forward Sequence:	TTGTGCTTGATCTAACCATG
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000888
Synonyms:	MI0000957; MIMAT0000888; MIR218a-2 rno-mir-218a-2; rno-miR-218a-5p
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)