

## Product datasheet for **RP300149**

### **MIR293 Rat qPCR Primer Pair (MI0012609)**

#### **Product data:**

|                   |                                                                                                                                                                                                                                                                                                                                                   |
|-------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                                 |
| Locus ID:         | 100314136                                                                                                                                                                                                                                                                                                                                         |
| Forward Sequence: | ACTCAAAGTGTGACACT                                                                                                                                                                                                                                                                                                                                 |
| Reverse Sequence: | GAACATGTCTGCGTATCTC                                                                                                                                                                                                                                                                                                                               |
| ACCN:             | MIMAT0012847                                                                                                                                                                                                                                                                                                                                      |
| Synonyms:         | MI0012609; MIMAT0012847; MIR293 rno-mir-293; rno-miR-293-5p                                                                                                                                                                                                                                                                                       |
| Components:       | 1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)                                                                                                                                                                                                                                                   |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec. |
| Storage:          | The primer mix is stable for one year from date of shipping. Store at -20 °C.                                                                                                                                                                                                                                                                     |



[View online »](#)