

Product datasheet for **RP300204**

MIR344b-1 Rat qPCR Primer Pair (MI0015472)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100526625
Forward Sequence:	AGTCAGGCTGCTGGTTAT
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0017811
Synonyms:	MI0015472; MIMAT0017811; MIR344b-1 rno-mir-344b-1; rno-miR-344b-5p
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20°C.



[View online »](#)