

## Product datasheet for **RP300309**

### **MIR466b-1 Rat qPCR Primer Pair (MI0006112)**

#### **Product data:**

|                   |                                                                                                                                                                                                                                                                                                                                               |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:     | qPCR Primer Pairs                                                                                                                                                                                                                                                                                                                             |
| Locus ID:         | 100314093                                                                                                                                                                                                                                                                                                                                     |
| Forward Sequence: | ATGTGTGTGTGTATGTCCA                                                                                                                                                                                                                                                                                                                           |
| Reverse Sequence: | GAACATGTCTGCGTATCTC                                                                                                                                                                                                                                                                                                                           |
| ACCN:             | MIMAT0005278                                                                                                                                                                                                                                                                                                                                  |
| Synonyms:         | MI0006112; MIMAT0005278; MIR466b-1 rno-mir-466b-1; rno-miR-466b-5p                                                                                                                                                                                                                                                                            |
| Components:       | 1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)                                                                                                                                                                                                                                               |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | The primer mix is stable for one year from date of shipping. Store at -20°C.                                                                                                                                                                                                                                                                  |



[View online »](#)