

## Product datasheet for **RP300374**

### MIR664-2 Rat qPCR Primer Pair (MI0003723)

#### Product data:

|                   |  |
|-------------------|--|
| Product Type:     | qPCR Primer Pairs  |
| Locus ID:         | 100314092  |
| Forward Sequence: | TGGCTGGGGAAAATGATTG  |
| Reverse Sequence: | GAACATGTCTGCGTATCTC  |
| ACCN:             | MIMAT0017229   |
| Synonyms:         | MI0003723; MIMAT0017229; MIR664-2 rno-mir-664-2; rno-miR-664-2-5p  |
| Components:       | 1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)  |
| Quality Control:  | The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak:95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60°C for 1 min; Stage 4: Melting curve: 95°C for 15 sec, 60°C for 15 sec, 95°C for 15 sec. |
| Storage:          | The primer mix is stable for one year from date of shipping. Store at -20°C.   |



[View online »](#)