

Product datasheet for **RP300400**

MIR7a-1 Rat qPCR Primer Pair (MI0000641)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100314148
Forward Sequence:	TGGAAGACTAGTGATTTTGT
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0000606
Synonyms:	MI0000641; MIMAT0000606; MIR7a-1 rno-mir-7a-1; rno-miR-7a-5p
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20 °C.



[View online »](#)