

Product datasheet for **RP300427**

MIR9b-1 Rat qPCR Primer Pair (MI0015430)

Product data:

Product Type:	qPCR Primer Pairs
Locus ID:	100526581
Forward Sequence:	TTCGGTTATCTAGCTTTAT
Reverse Sequence:	GAACATGTCTGCGTATCTC
ACCN:	MIMAT0017835
Synonyms:	MI0015430; MIMAT0017835; MIR9b-1 rno-mir-9b-1; rno-miR-9b-5p
Components:	1 vial of lyophilized qSTAR miRNA primer mix (1 nmol each primer, sufficient for 200 reactions)
Quality Control:	The primer mix has been tested to generate satisfactory qPCR data on ABI 7900HT by using the following PCR program: Stage 1: Activation: 50 °C for 2 min; Stage 2: pre-soak: 95 °C for 10 min; Stage 3: Denaturation: 95 °C for 15 sec, Annealing: 60 °C for 1 min; Stage 4: Melting curve: 95 °C for 15 sec, 60 °C for 15 sec, 95 °C for 15 sec.
Storage:	The primer mix is stable for one year from date of shipping. Store at -20 °C.

[View online »](#)