

Product datasheet for **RR200007**

Sinhcaf (NM_001134711) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
 Product Name: Sinhcaf (NM_001134711) Rat Tagged ORF Clone
 Tag: Myc-DDK
 Symbol: Sinhcaf
 Synonyms: Fam60a; MGC187734
 Vector: pCMV6-Entry (PS100001)
 E. coli Selection: Kanamycin (25 ug/mL)
 Cell Selection: Neomycin
 ORF Nucleotide Sequence: >RR200007 representing NM_001134711
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTTTGGTTTTACAAGCCAAAGATGTACCGAAGTATAGAAGGCTGCTGTATCTGCAGAGCCAAGTCTT
 CCAGCTCTCGGTTACGGACAGCAAGCGCTATGAGAAGGACTTCCAGAGCTGCTTTGGGCTGCATGAGAC
 TCGCTCAGGAGACATCTGCAACGCCTGTGTGCTGCTTGTGAAAAGATGGAAGAAGCTGCCGCAGGATCC
 AAAAAAACTGGAATCAGTCGTAGATGCAAGAGCAGGCCCCAGTCTAAAGACAACACTGAAACCAAGA
 AAGTGAAGACTCTGTCTGGAAACAGGATGAAAAGCAGCCAGATCAGTAAGCTGCAGAAGGAGTTAAACG
 CCACAATTCTGATGCTCATAGCACCCACTCAAGTGCTCTCCAGCCCAGTCTCCATGCTATAGTAACCG
 TCAGACGAGGGCTCAGATACAGAGATGGCTTCCAGCTCTAATAGAACTCCGGTTTTTCTTCTTAGATC
 TTACCTACTGGAAAAGACAGAAAATATGTTGTGGGATCATCTATAAGGGCCGTTTTGGGGAAGTCTCAT
 CGACACTCATCTTTCAAGCCTTGCTGCAGCAGTAAGAAAAGCTGCTGCCGAGAAGCCAGAGGAGCAGGCC
 CCAGCGCCTTGCCCATCTCCACTCAGGAGTGG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200007 representing NM_001134711
 Red=Cloning site Green=Tags(s)

MFGFHKPKMYRSIEGCCICRAKSSSSSRFTDSKRYEKDFQSCFGLHETRSGDICNACVLLVKRWKLPAGS
 KKNWNHVVDARAGPSLKTTLKPKVKTL SGNRMKSSQISKLQKEFKRHNSDAHSTTSSASPAQSPCYSNQ
 SDEGSDTEMASSNRTPVF SFLDLTYWKRQKICCGI IYKGRFGEVLIDHLFKPCCSSKAAAEEKPEEQA
 PAPLP ISTQEW

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI


Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN: NM_001134711

ORF Size: 663 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001134711.1 , NP_001128183.1
RefSeq Size:	2781 bp
RefSeq ORF:	666 bp
Locus ID:	686611
Cytogenetics:	4q44
MW:	24.8 kDa