

Product datasheet for **RR200011**

Zfp385a (NM_001135089) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Zfp385a (NM_001135089) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Zfp385a
Synonyms: Zfp385
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200011 representing NM_001135089
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCAGCCCCGATGGACCTGAAGCAGATCCTGCCCTTCCCACTGGAGCCAGCCCCAACCTGGGCCTCT
TCAGAACTACAGACAATGGACCCTGTACAGAAGGCTGTGCTCTCCACACATTTGGAGGACCCTTGCT
CAAGACCAAGCGCCAGTCATTTCTGTAAATGTCTGTCAAATCCGCTTCAATTCTCAGAGCCAGGCTGAG
GCGCACTACAAGGTAATCGCCATGCCGAAGAGTCAAAGGCATCGAAGCTGCCAAAACCCGAGGCAGGG
AGCCTAGTGTCCGGGAATCAGGAGATCCGGCTCCAGCAGGCAGCACCCCTCCAAGTGGGGATGGCGTAGC
CCCTCGTCCAGTTTCCATGGAGAATGGCCTGGGTCCAGCCCCAGGATCCCCAGAGAAACAGCCGGGCTCC
CCATCTCCTCCAGTGTTCAGAGTCTGGACAGGGTGTAAACCAAGGTGAAGGGGGGACTTCAGTCCCAG
CTTCCCTGCCTGGGGGTAGCAAGGAAGAAGAGGAGAAAGTAAGCGTCTGCTCTACTGTGCACTGTGCAA
GGTGGCTGTGAATCCCTGTCCCAGCTTGAGGCACATAACAAAGGTAAGCACAAGACAATTTGGAG
GCCCCAAGTGGGCTGGGACCCATCAAAGCTTACCCTCGTTGGGGCTCCCACTCCTGGGAACCCAGAGG
CTCCTGCCAGGACCGAACCTTCCACTGTGAGATCTGCAATGTCAAGGTCAATTCGGAGTCCAGCTGAA
ACAGCACATCTCCAGCAGGAGGCCACCGAGATGGCGTGGCTGGGAAGCCCAACCCTCTACTGAGCCGGCAC
AAGAAACCTAGGGGCGCTGCAGAGCTGGCGGGCAGCTGACTTTCTCCAAGGAGCTGCCAAGTCCCTGG
CCGGTGGCTGCTCCCCAGCCCCCTAGCGGTGGCTGCGGTGATGGCCGCTGCAGCAGGCTCTCCGCTGTC
CCTGCGCCAGCTCCAGCTGCACCTTTCTGCAGGGACCACCGATCACACACCCTCTACTCCACCCTGCC
CCAGGACCCATCCGAATGCGCACGGACCCATCCTTTCTCCCCCTAC

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200011 representing NM_001135089
 Red=Cloning site Green=Tags(s)

MQPPMDLKQILPFPLEPAPTLGLFSNYSTMDPVQKAVLSHTFGGPLLKTKRPVISCNVCQIRFNSQSQAE
 AHYKGNRHARRVKGIEAAKTRGREPSVRESGDPAPAGSTPPSGDGVAPRPVSMENGLGPAPGSPEKQPGS
 PSPPSVPESGQGVTKGEGGTSVPASLPGGSKEEEEKAKRLLYCALCKVAVNSLSQLEAHNKGTKHKHTILE
 ARSGLGPIKAYPRLGPPTPGEPEAPAQDRTFHCEICNVKVNSEVQLKQHISSRRHRDGVAGKPNPLLSRH
 KKPRGAAELAGTLTFSKELPKSLAGLLPSLAVAAMAAAAGSPLSLRPAPAAPLLQGPPITHPLLHPA
 PGPVRTAHGPILFSPY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001135089

ORF Size:	1098 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001135089.1 , NP_001128561.1
RefSeq Size:	2344 bp
RefSeq ORF:	1101 bp
Locus ID:	685474
Cytogenetics:	7q36
MW:	38.3 kDa