

## Product datasheet for **RR200025**

### Wtap (NM\_001113542) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Wtap (NM_001113542) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Wtap
Synonyms:	RGD1563824
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR200025 representing NM_001113542. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGACCAACGAAGAACCTCTTCTAAAAAGGTCCGACTGAGTGAACAGACTTCAAAGTTATGGCACGG
GATGAGTTAATTCTGAGATGGAAACAATATGAAGCATATGTTCAAGCTTTGGAGGGAAAGTACACAGAC
CTTAATCAAACGATGTGACTGGCTTAAGGGAATCTGAAGAAAACTAAAGCAGCAACAGCAGGAATCT
GCACGCAGGGAGAACATCCTTGTATGCGGCTAGCAACCAAGGAACAGGAGATGCAAGAGTGCACCACT
CAAATCCAGTACCTCAAGCAAGTTCAGCAGCCGAGTGTAGCCCAACTGAGATCAACAATGGTAGACCCA
GCGATCAACTTGTCTTCTAAAAATGAAAGGTGAACTGGAACAGACTAAAGACAACTGGAACAAGCC
CAAAATGAACTGAGTGCCTGGAAGTTTACGCCTGATAGCCAAACAGGCAAAAAATTAATGGCGAAGTGT
CGAATGCTTATCCAGGAGAACCAAGAGCTTGAAGGCAGCTGTCCAGGGCCGTATTGCACAGCTTGAA
GCCGAGTTGGCTTTACAGAAGAAATACAGTGAGGAGCTTAAAAGCAGTCAGGATGAACTGAATGACTTC
ATCATTCAACTTGACGAAGAAGTAGAGGGTATGCAGAGCACCATTCTAGTTCTTCAGCAACAGTTGAAG
GAGACTCGACAGCAGTTGGCACAGTACCAGCAGCAGTCTCAAGCTTCTGCTCAAGTACCAGCAGG
ACTACATCTTCTGAACCTGTGGATCAGGCAGAGGCCACAAGCAAGACTGCAGTCGTCTGGCAAATGGA
CCAAGTAATGGCAGTCTCTCCCGCCAGAGGACGTCTGGGTCTGGATTTCACAGGGAGGGGAGCACTCCT
GAGGATGACTTTCTTCTTCTCAGGAAACGGTAATAAGGCCTCCAACAGCTCAGAGGAGAGAAGTGGC
AGAGGTAGTAGTTACATAAAACCACTCAGTGCAGGGTATGAAAGTGTAGACTCTCCACGGGCAGTGAA
AACTCTCTCACACATCACTCAAATGACACAGACTCCAGTCATGACCCTCAAGAGGAGAAAGCAGTGAGT
GGGAAAGGTAATCGAACTGTGGTTCCCGCCATGTTTCAAGATGGCTTGGACTCAAGTGTAAATGTACAG
GGTTCAATTTTG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```




[View online >](#)

**Protein Sequence:** >Peptide sequence encoded by RR200025  
 Blue=ORF Red=Cloning site Green=Tag(s)


MTNEEPLPKKVRLSEDFKVMARDELILRWKQYEAYVQALEGKYTDLNSNDVTGLRESEELKQQQQES  
 ARRENILVMRLATKEQEMQECTTQIQYLKQVQQPSVAQLRSTMVDPAINLFFLKMKGELEQTKDKLEQA  
 QNELSAWKFTPDSQTGKKLMAKCRMLIQENQELGRQLSQGRIAQLEAELALQKKYSEELKSSQDELNDF  
 I IQLDEEVEGMQSTILVLQQQLKETRQQLAQYQQQQSQASAPSTRRTTSSEPVQAEATSKDCSRLANG  
 PSNGSSSRQRTSGSGFHREGSTPEDDFPSSSGNGNKASNSSEERTGRGSSYINPLSAGYESVDSPTGSE  
 NSLTHHSNDTSSHDPQEEKAVSGKGNRTVGSRHVQNGLDSSVNVQGSIL  
 TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001113542

**ORF Size:** 1185 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001113542.1</a> , <a href="#">NP_001107014.1</a>
<b>RefSeq Size:</b>	2028 bp
<b>RefSeq ORF:</b>	1188 bp
<b>Locus ID:</b>	499020
<b>Cytogenetics:</b>	1q11
<b>MW:</b>	44.1 kDa