

## Product datasheet for RR200037L4

### Dnmt3a (NM\_001003958) Rat Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Dnmt3a (NM_001003958) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Dnmt3a
Synonyms:	MGC109195
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200037).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.



[View online »](#)



Locus ID: 444984  
UniProt ID: [Q1LZ53](#)  
Cytogenetics: 6q14  
Gene Summary: member of a family of DNA methyltransferases [RGD, Feb 2006]