

Product datasheet for **RR200054**

Svbp (NM_001038992) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Svbp (NM_001038992) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Svbp
Synonyms: Ccdc23; RGD1311232
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200054 representing NM_001038992
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGATCCACCTGCCCGAAAGAAAATCCAAAGTTAAAGAACCTGCCTTCAGAGTGGAGAAGGCTAAAC
AGAAATCTGCCAGCAGGAGCTGAAGCAAAGACAAGAGCAGAGATCTATGCTCTCAACAGATCATGAC
GGAGCTGGAGCAGCAGCAGTTTGATGAGTTCTGTAAAGCAGATGCAGCCCCCGGGGAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR200054 representing NM_001038992
Red=Cloning site Green=Tags(s)

MDPPARKEKSKVKEPAFRVEKAKQKSAQQELKQRQRAEIYALNRVMTELEQQQFDEFCKQMPPGE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-Mlul




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001038992
 ORF Size: 198 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001038992.1 , NP_001034081.1
RefSeq Size:	688 bp
RefSeq ORF:	201 bp
Locus ID:	362578
UniProt ID:	Q4KLG3
Cytogenetics:	5q36
MW:	7.8 kDa
Gene Summary:	Enhances the tyrosine carboxypeptidase activity of VASH1 and VASH2, thereby promoting the removal of the C-terminal tyrosine residue of alpha-tubulin. Also required to enhance the solubility and secretion of VASH1 and VASH2.[UniProtKB/Swiss-Prot Function]