

## Product datasheet for **RR200077**

### Tbc1d14 (NM\_001012152) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tbc1d14 (NM\_001012152) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Tbc1d14  
**Synonyms:** SRF-2; UR-NR#2  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR200077 representing NM\_001012152  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGAATATGAAGACAAGGCTGGAAGACCAAGCAGGCCGCCCTCTCAAAGCAGAATGTGAGGAAGAACC  
 TGGATTTT**GAGCCACTGTCCACCACAGCGCTCATCCTTGAGGACAGACCCGCGAATCTCCAGCAAAGCC**  
**AGCTGAAGAAGCTCAAAGCACAGACAGCAGTATGAAGAGATGGTGGTGCAGGCTAAGAAACGAGAGCTT**  
**AAAGAAGCCAGCGCGGAAGAAACAGCTGGAGGAAAGGTGCAAGTTGAGAGAGCATTGAAACGCCG**  
**TCCTCACCTGGAATAATGAGATCTTGCCTAACTGGGAAACCATGTGGTGCCTAAGAAGTTTCGAGACTT**  
**GTGGTGGCAGGGAATCCCTCCGAGTGTAAAGAGCAAGTTTGGAGCTTGGCCATTGGCAACGAGCTGAAC**  
**ATCACTCACGAGCTCTTTGACATCTGTCTTGTCTCGGGCAAGGAGCGGTGGCGTCCCTCAGCACTGGAG**  
**GCTCAGAAGTGGAGAATGAAGATGCTGGGTTTTCAGCAGCAGACAGAGAGGCCAGTCTGGAGCTTATTA**  
**ACTGGACATTTCTAGAACATTTCTAATCTCTGCATTTCCAGCAAGTGGTCCATATCATGACATGTTG**  
**CACAGTATTTGGGCGCTTATACTTGTACCAGCGGATGTGGTTATGTT**CAGGGCATGCTCTTTATAG****  
**CTGCAGTGTGATCTTGAACCTCGATACTGCAGATGCCTTATTGCTTTTTCTAATCTTTGAATAAACC**  
**CTGTCAAATGGCGTTTTTCCGAGTGGACCATGGCCTTATGTTGACTTACTTTGCTGCATTTGAGGTATTC**  
**TTTGAAGAAAACCTTGCCGAAATTGTTTGCACATTTCAAGAAAAACAACCTAACTGCAGACATCTACCTAA**  
**TTGATTGGATTTTACCTTGACAGTAAGTCCCTGCCCTTGGACTTGGCCTGCCGATCTGGGATGTATT**  
**CTGTCGAGATGGAGAGGAATTTCTCTCCGACAGCTCTGGGCATCCTGAAGCTTTTGGAGATACCTG**  
**ACCAGGATGGACTTCATCCACAGTGCCAGTTTCTGACCAGGCTGCCTGAGGACCTGCCAGCAGACGACG**  
**TGTTCCGACGATCTCTACAGTCCAGATGCAAAGCCGCAACAAGAAATGGGCACAGGTGCTGTCTGCACT**  
**GCAGAAAGACAGCCGGAAATGGAGAAAGGAAGCCCATCACTCAGACAC**

**ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT**  
**ACAAGGATGACGACGATAAGGTTTAA**



[View online »](#)

**Protein Sequence:** >RR200077 representing NM\_001012152  
 Red=Cloning site Green=Tags(s)

MEYEDKAGRPSRPPSPKQNVKLNDFEPLSTTALILEDSPANLPAKPAAEAQKHRQQYEEMVQAKKREL  
 KEAQRKKQLEERCKVEESIGNAVLTWNNEILPNWETMWCSSKVRDLWWQGIPPSVRGKVWSLAIGNELN  
 ITHLFDICLARAKERWRSLSSTGGSEVENEDAGFSAADREASLELIKLDISRTFPNLCIFQQGGPYHDM  
 HSILGAYTCYRPDVGIVQGMFSIAAVLILNLDTADAFIAFSNLLNKPCQMAFFRVDHGLMLTYFAAFEV  
 FEENLPKLFHFHFKNNLTADIYLDWIFTLYSKSLPLDLACRIWDVFCRDGEEFLFRALTALGILKLFEDIL  
 TRMDFIHSAQFLTRLPELDPADDVFAAISTVQMQRNKKWAQVLSALQKDSREMEKGSPLRHH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001012152

<b>ORF Size:</b>	1239 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001012152.2</a> , <a href="#">NP_001012152.2</a>
<b>RefSeq Size:</b>	3859 bp
<b>RefSeq ORF:</b>	1242 bp
<b>Locus ID:</b>	360956
<b>UniProt ID:</b>	<a href="#">Q5CD77</a>
<b>Cytogenetics:</b>	14q21
<b>MW:</b>	47.6 kDa
<b>Gene Summary:</b>	Negative regulator of starvation-induced autophagosome formation.[UniProtKB/Swiss-Prot Function]