

Product datasheet for **RR200078L4V**

Tbc1d14 (NM_001113365) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tbc1d14 (NM_001113365) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tbc1d14
Synonyms:	SRF-2; UR-NR#2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001113365
ORF Size:	2142 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200078).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001113365.1 , NP_001106836.1
RefSeq Size:	4948 bp
RefSeq ORF:	2145 bp
Locus ID:	360956
UniProt ID:	Q5CD77
Cytogenetics:	14q21
Gene Summary:	Negative regulator of starvation-induced autophagosome formation.[UniProtKB/Swiss-Prot Function]



[View online »](#)