

Product datasheet for **RR200079**

Dnaja3 (NM_001038595) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Dnaja3 (NM_001038595) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dnaja3
Synonyms:	Tid-1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>NM_001038595 ORF sequence, RR200079 may differ due to SNPs. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCTGCGCGGTGTTCCCCACGCTGGTTAGGGTGGCGGTGGGGACCCCGGCTGCCGGCTGCAGCG
GGGAGAGGGGTCCAGCAGCCCCAGGGGGCGTAGTGGCAGCATCCCTGTGTCGAAAGCTGTGTGTCTCT
GCCTTCGGGCTTCCATGGGCGCTCACGCCCCCGGGCGCTGCTGACATTGAGACCCGGTGTCCGCCTC
ACAGGAACAAAAAGTTTCCCTTTTGTGTTACTGCTTCTTCCACACAAGTGCTTCTTTGGCCAAAGAT
GATTATTACCAGATACTGGGAGTACCCGAAATGCCAGCCAGAAGGATATCAAGAAAGCCTATTATCAG
CTTGCCAAAGAAATATCATCCAGACACAAATAAGGATGATCCCAAAGCCAAGGAGAAGTTTTCCAGCTG
GCAGAAGCTTATGAGGTGTTGAGTGACGAGGTGAAGAGGAAGCAGTATGATGCTTACGGCTCCGCTGGC
TTTGACCTGGAGCCAGCAGCTCTGGGCAGGGCTACTGGAGAGGAGGTCCTTCTGTTGACCCGAGGAG
CTATTCAGGAAGATCTTTGGAGAGTTCTCATCTTCTCCTTTTGGTGATTCCAGAATGTGTTTGATCAG
CCTCAGGAGTACATCATGGAGTTGACATTCAATCAAGCCGCCAAGGGTGTCAACAAGAGTTCACCGTG
AACATCATGGACACCTGTGAGCGCTCGATGGCAAGGGGAACGAGCCTGGAACCAAGTGCAGCATTGT
CACTACTGCAGCGGCTCAGGCATGGAAACCATCAATACAGGGCCTTTTGTGATGCGTTCCACGTTCGG
AGATGTGGTGGACGGGGCTCCATCATCACAAACCCTTGTGGTCTGCAGAGGAGCAGGACAAGCCAAG
CAGAAGAAGCGAGTGACAGTCCCTGTGCCTGCAGGAGTTGAAGATGGCCAGACTGTAAGGATGCCTGTG
GGAAAACGAGAAATTTTTGTACATTACAGGTGCAGAAGAGCCCTGTGTTCCGGAGGGACGGCGCAGAC
ATCCACTCCGACCTCTTTATTTCAATAGCTCAAGCTATTCTTGGGGGTACAGCCAAAGCGCAGGGCCTG
TATGAGACAATCAATGTGACGATCCCTGCTGGGATTACAGACAGACCAGAAGATTCGCTTGACTGGGAAA
GGCATCCCCGTATTAATAGCTATGGCTATGGTGACCACTACATTACATTAAGATCAGAGTTCCAAAG
AGGCTAAGCAGTAGGCAGCAGAACCTGATCCTGAGCTATGCTGAGGATGAGACCGATGTGGAGGGGACA
GTGAATGGAGTACCCACACAAGCACTGGTGGCAGAACCATGGATAGCTCCACAGGAAGCAAGGATAGG
CGAGAGGCTGGGAGGACAATGAGGGATTCTTTCCAAACTTAAGAAAATATTTACTTCC
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```




Protein Sequence: >Peptide sequence encoded by RR200079
 Blue=ORF Red=Cloning site Green=Tag(s)


MAARCSPRWFRVAVGTPRLPAAAGRGVQQPQGGVVAASLCKKLCVSAFGLSMGAHGPRALLLRPGVRL
 TGTKSFPFVCTASFHTSASLAKDDYYQILGVPRNASQKDIKKAYQLAKKYHPDTNKDDPKAKEKFSQL
 AEAYEVLSDVEVSRKQYDAYGSAGFDPGASSGQGYWRGGPSVDPEELFRKIFGEFSSSPFGDFQNVFDQ
 PQEYIMELTFNQAAKGVNKEFTVNMIMDTCERCDCGKGNPVGTKVQHCHYCSGSGMETINTGPFVMRSTCR
 RCGGRGSIITNPCVFCRGAGQAKQKRVTVVPVAGVEDGQTVRMPVKGREIFVTFRVQKSPVFRRDGAD
 IHSDLFISIAQAILGGTAKAQGLYETINVTIPAGIQTDQKIRLTGKGIPRINSYGYGDHYIHIKIRVPK
 RLSSRQQLILSYAEDETDVEGTVNGVTHSTGGRTMDSSTGSKDRREAGEDNEGFLSKLKKIFTS
 TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001038595

ORF Size:	1440 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001038595.2
RefSeq Size:	2645 bp
RefSeq ORF:	1443 bp
Locus ID:	360481
Cytogenetics:	10q12
MW:	52.4 kDa