

Product datasheet for RR200083L4

Gmcl1 (NM_001010956) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gmcl1 (NM_001010956) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Gmcl1
Synonyms:	Gcl
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200083).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_001010956
ORF Size:	1479 bp



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001010956.3 , NP_001010956.1
RefSeq Size:	1952 bp
RefSeq ORF:	1482 bp
Locus ID:	312516
Cytogenetics:	4q34
Gene Summary:	This gene encodes a nuclear envelope protein that contains a broad-complex, tramtrack and bric a brac/ Pox virus and Zinc finger (BTB/POZ) protein-protein interaction domain. In mouse, this gene is expressed at high levels in the testis, and null mutant mice exhibit morphologically abnormal sperm that have chromatin condensation defects. These defects result in significantly impaired fertility, demonstrating a role for this gene in male germ cell development in mouse. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2015]

Product images:



Circular map for RR200083L4