

Product datasheet for **RR200085L4V**

Cdc14a (NM_001107718) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdc14a (NM_001107718) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cdc14a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001107718
ORF Size:	1878 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200085).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001107718.2 , NP_001101188.2
RefSeq Size:	2613 bp
RefSeq ORF:	1881 bp
Locus ID:	310806
Cytogenetics:	2q41



[View online »](#)

Gene Summary:

The protein encoded by this gene is a dual-specificity phosphatase that preferentially dephosphorylates cyclin dependent kinase substrates to regulate the cell cycle. In human cell lines, this protein localizes to interphase chromosomes, and depletion of the transcript results in centrosome separation and cytokinesis defects. In mouse, the protein localizes to the nucleus of prophase I arrested oocytes and then becomes dispersed in meiotically competent oocytes. Knockdown of the protein delays exit from metaphase I and results in eggs with chromosomal abnormalities and elevated aneuploidy, demonstrating a function in regulation of meiosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2015]