

Product datasheet for **RR200104**

Pbx1 (NM_001134862) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pbx1 (NM_001134862) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pbx1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR200104 representing NM_001134862 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGACGAGCAGCCGAGGCTGATGCATTCCCACGCTGGGGTCGGGATGGCCGGACACCCCGGCTGTCCC
AGCACTTGCAGGATGGGCGGAGGGACCGAGGGGAGGGCGGGAGGAAGCAGGACATCGGGACATTTT
ACAGCAAATTATGACCATCACAGACCAGAGTTTGGATGAAGCGCAGGCCAGAAAACATGCTTTAACTGC
CACAGAATGAAGCCTGCCTTGTTAATGTGTTGTGTGAAATCAAAGAAAAACAGTTTTGAGTATTCGGG
GAGCCCAAGAAGAGGAGCCACAGACCCACAGCTCATGCGACTGGACAACATGCTGCTAGCAGAAGGGT
GGCGGGGCTGAGAAGGGCGGAGGCTCGGCAGCGGCAGCGCAGCGGAGCTTCTGGGGTGCAGGT
TCAGACAACCTCAGTGGAGCATTCCGATTACAGAGCCAACTCTCACAGATCAGACAAATCTACCACAG
AGCTGGAGAAGTATGAACAGGCATGCAATGAATTCACCACCCACGTGATGAACCTCCTTCGGGAGCAAAG
CCGGACCAGACCCATCTCTCCGAAGGAGATCGAGCGGATGGTCAGCATTATTCACCGCAAGTTCAGCTCC
ATCCAGATGCAGCTGAAGCAAAGCACGTGCGAGGCAGTCATGATCCTGCGCTCCCGTTTCTGGATGCGA
GGCGGAAGAGACGGAATTTCAACAAGCAAGCCACAGAAATCTGAATGAATATTTCTATCCCATCTCAG
CAACCTTACCCAGTGAGGAAGCCAAAGAGGAGTTAGCCAAGAAGTGGCGCATCACAGTCTCCAGGTA
TCAAACCTGGTTTGAAATAAGCGAATCCGGTACAAGAAGAACATAGGTAAATTTCAAGAGGAAGCCAATA
TTTATGCTGCCAAAACGGCTGTCACTGCCACCAATGTGTGAGCCCATGGAAGCCAAGCTAACTCGCCCTC
CACTCCAACTCAGCGGGTTCTTCCAGTCTTTTAAACATGTCAAACCTCTGGAGATTTGTTTCATGAGCGTG
CAGTCACTCAATGGGATTTTACCAAGGGGCCAGGTTGGAGCCAACGTGCAATCACAGGTGGATACCC
TTCGCCATGTTATCAGCCAGACAGGAGGATACAGTGACGGACTCGCAGCCAGTCAGATGTACAGTCCGCA
GGGCATCAGTGCTAATGGAGTTGGCAGGATGCTACTACCCCTTCATCAGTGACCTCCCTACAGAAGGC
CCTGGCAGTGTCACTCTGATACCTCCAAC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200104 representing NM_001134862
 Red=Cloning site Green=Tags(s)


MDEQPRLMHSAGVGMAGHPGLSQHLQDGAGGTEGEGGRKQDIGDILQQIMTITDQSLDEAQARKHALNC
 HRMKPALFNVLCEIKEKTVLSIRGAQEEPTDPQLMRLDNMLLAEGVAGPEKGGGSAASAAAAAASGGAG
 SDNSVEHSDYRAKLSQIRQIYHTELEKEYEQACNEFTTHVMNLLREQSRTRPISPKEIERMVSIIHRKFSS
 IQMQLKQSTCEAVMILRSRFLDARRKRRNFNKQATEILNEYFYSHLSNPYPSEEAKEELAKKCGITVSQV
 SNWFGNKRIRYKKNIGKFQEEANIYAAKTAVTATNVSAHGSQANSPSTPNASGSSSSFNMSNSGDLFMSV
 QSLNGDSYQGAQVGANVQSQVDTLRHVISQTGGYSDGLAASQMYSPQGISANGGWQDATTSSVTSPTTEG
 PGSVHSDTSN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:	NM_001134862
ORF Size:	1290 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001134862.2
RefSeq Size:	2896 bp
RefSeq ORF:	1293 bp
Locus ID:	304947
Cytogenetics:	13q24
MW:	46.6 kDa