

Product datasheet for **RR200105**

Pbx1 (NM_001100681) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pbx1 (NM_001100681) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pbx1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR200105 representing NM_001100681 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACGAGCAGCCGAGGCTGATGCATTCCCACGCTGGGGTCGGGATGGCCGGACACCCCGCCTGTCCC
AGCACTTGCAGGATGGGCGGGAGGGACCGAGGGGGAGGGCGGGAGGAAGCAGGACATCGGGGACATTTT
ACAGCAAATTATGACCATCACAGACCAGAGTTTGGATGAAGCGCAGGCCAGAAAACATGCTTTAACTGC
CACAGAATGAAGCCTGCCTTGTTAATGTGTTGTGAAATCAAAGAAAAACAGTTTTGAGTATTCGGG
GAGCCCAAGAAGAGGAGCCACAGACCCACAGCTCATGCGACTGGACAACATGCTGCTAGCAGAAGGGT
GGCGGGGCTGAGAAGGGCGGAGGCTCGGCAGCGGCAGCGCAGCAGCGCAGCTTCTGGGGTGCAGGT
TCAGACAACCTCAGTGGAGCATTCCGATTACAGAGCCAACTCTCACAGATCAGACAAATCTACCACACAG
AGCTGGAGAAGTATGAACAGGCATGCAATGAATTCACCACCCACGTGATGAACCTCCTTCGGGAGCAAAG
CCGGACCAGACCCATCTCTCCGAAGGAGATCGAGCGGATGGTCAGCATTATTCACCGCAAGTTCAGCTCC
ATCCAGATGCAGCTGAAGCAAAGCACGTGCGAGGCAGTCATGATCCTGCGCTCCCGTTTCTGGATGCGA
GGCGGAAGAGACGAATTTCAACAAGCAAGCCACAGAAATCTGAATGAATTTTCTATTCCCATCTCAG
CAACCTTACCCAGTGAGGAAGCCAAAGAGGAGTTAGCCAAGAAGTGGCGCATCACAGTCTCCAGGTA
TCAAACCTGGTTTGGAAATAAGCGAATCCGGTACAAGAAGAACATAGGTAATTTCAAGAGGAAGCCAATA
TTTATGCTGCCAAAACGGCTGTCACTGCCACCAATGTGTGAGCCCATGGAAGCCAAGCTAACTCGCCCTC
CACTCCAACTCAGCGGGTGGATACCCTTCGCCATGTTATCAGCCAGACAGGAGGATACAG

ACGCGTACGCGGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200105 representing NM_001100681
 Red=Cloning site Green=Tags(s)

MDEQPRLMHSAGVGMAGHPGLSQHLQDGAGGTEGEGGRKQDIGDILQQIMTITDQSLDEAQARKHALNC
 HRMKPALFNVLC EIKEKTVLSIRGAQEEPTDPLMRLDNMLLAEGVAGPEKGGSSAAAAAASGGAG
 SDNSVEHSDYRAKLSQIRQIYHTELEKYEQACNEFTTHVMNLLREQSRTRPISPKEIERMVSIIHRKFSS
 IQMQLKQSTCEAVMILRSRFLDARRKRRNFNKQATEILNEYFYSHLSNPYPSEEAKEELAKKCGITVSQV
 SNWFGNKRIRYKKNIGKFQEEANIYAAKTAVTATNVSAHGSQANSPSTPNSAGGYPSPCYQPDRIQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001100681

ORF Size:	1041 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001100681.2
RefSeq Size:	1869 bp
RefSeq ORF:	1044 bp
Locus ID:	304947
Cytogenetics:	13q24
MW:	38.4 kDa