

## Product datasheet for **RR200106**

### Galnt9 (NM\_001107151) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Galnt9 (NM\_001107151) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Galnt9  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR200106 representing NM\_001107151  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGATCGCC

ATGGAGGTACTCCCTGCTCCCGTGTGGCCACATTGAGCGTACCAAAAAGCCATATAACAATGACATCG  
ACTACTACGCCAAGCGTAACGCGCTGAGGGCCGCGAGGTGGATGGATGATTTTAAGTCCCACGTCTA  
CATGGCATGGAATATCCCATGACAACCCAGGAGTGGACTTTGGGGATGTGCTGAGAGGCTGGCCCTG  
CGTCAGAAAGTTAAAATGCCGTAGCTTCAAATGGTACCTGGAGAACGTGTACCCCGAGATGAGGATCTACA  
ACAACACCCTCACGTACGGAGAGGTGAGAAACAGCAAAGCCAGCGGCTACTGCTTGGACCAGGGGCGAGA  
AGACGATGACCGAGCCATTCTTTACCCTTGCCATGGCATGTCCTCACAGCTGGTGCCTACAGTGTGAT  
GGCCTCTGCAGCTGGGCCCCCTGGGCTCTACCGCTTTCTGCCAGACTCCAAGTGCCTGGTAGATGACG  
GCCGGGGACGCATGCCACCCTAAAGAAGTGTGAGGATGTAGCCAGGCCTACGCAGAGACTGTGGGACTT  
CACCCAGAGTGGCCCCATTGTGAGCCGAGACACAGGTCGGTGTCTGGAGGTAGAGATGTCTAAAGATGCC  
AACTTTGGGCTGAGGTTGGTGGTGCAGCGCTGCTCAGGGCAGAAGTGGATGATCCGAAACTGGATCAAAC  
ACGCCAAGCCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR200106 representing NM\_001107151  
Red=Cloning site Green=Tags(s)

MEVLPCSRVAHIERTKKPYNNDIDYYAKRNALRAAEVWMDDFKSHVYMAWNI<sup>1</sup>PMTNPGVDFG<sup>2</sup>VDV<sup>3</sup>SERLAL<sup>4</sup>  
 RQ<sup>5</sup>LKCR<sup>6</sup>SFKW<sup>7</sup>YLENV<sup>8</sup>PEM<sup>9</sup>RIY<sup>10</sup>NN<sup>11</sup>TYGE<sup>12</sup>VR<sup>13</sup>NSK<sup>14</sup>AS<sup>15</sup>GY<sup>16</sup>CLD<sup>17</sup>QGAED<sup>18</sup>DD<sup>19</sup>RA<sup>20</sup>ILY<sup>21</sup>PCH<sup>22</sup>GM<sup>23</sup>SS<sup>24</sup>QL<sup>25</sup>V<sup>26</sup>RY<sup>27</sup>SAD<sup>28</sup>  
 GLL<sup>29</sup>QL<sup>30</sup>GPL<sup>31</sup>G<sup>32</sup>STAF<sup>33</sup>LP<sup>34</sup>DSK<sup>35</sup>CL<sup>36</sup>VDD<sup>37</sup>GR<sup>38</sup>MP<sup>39</sup>TL<sup>40</sup>KK<sup>41</sup>CE<sup>42</sup>D<sup>43</sup>VAR<sup>44</sup>P<sup>45</sup>TQ<sup>46</sup>RL<sup>47</sup>W<sup>48</sup>DF<sup>49</sup>TQ<sup>50</sup>SG<sup>51</sup>PI<sup>52</sup>VR<sup>53</sup>SD<sup>54</sup>TGR<sup>55</sup>CLE<sup>56</sup>VEM<sup>57</sup>SK<sup>58</sup>DA<sup>59</sup>  
 NFGL<sup>60</sup>RL<sup>61</sup>V<sup>62</sup>V<sup>63</sup>QR<sup>64</sup>CS<sup>65</sup>G<sup>66</sup>Q<sup>67</sup>K<sup>68</sup>WM<sup>69</sup>IR<sup>70</sup>N<sup>71</sup>W<sup>72</sup>IK<sup>73</sup>H<sup>74</sup>AK<sup>75</sup>P


TRTRPLEQKLISEEDLANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM\_001107151

ORF Size: 711 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

|                               |   |
|-------------------------------|---|
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u>NM_001107151.2, NP_001100621.1</u>   |
| <b>RefSeq Size:</b>           | 1546 bp   |
| <b>RefSeq ORF:</b>            | 714 bp  |
| <b>Locus ID:</b>              | 304571  |
| <b>Cytogenetics:</b>          | 12q16   |
| <b>MW:</b>                    | 27.3 kDa  |