

Product datasheet for **RR200116**

Sp7 (NM_001037632) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sp7 (NM_001037632) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Sp7
Synonyms: Osx
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200116 representing NM_001037632
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGCTGACAGCAGCCTGCAGCAAGTTTGGTGGCTCCAGCCCTCTGCGAGACTCAACAGCCCTGGGAAAAG
GAGGCACAAAGAAGCCATACACTGACCTTTCAGCCCCAAAACCATGGGGGACGCCTACCGAGCTCCCTT
CTCAAGCACCAATGGACTCCTCTCCCTGCAGGCAGTCCCCGGCCCCGGCTCTGGCTATGCCAATGAC
TACCCACCCTTTCCCACTCATTTCCTGGGCCACTGGTGCCCAAGACCCGGGGCTGCTGGTGCCGAAGG
GGCACAGCTCTTCTGACTGCCTGCCTAGTGTCTACACGTCCTTGGATATGTCCCATCCCTACGGCTCCTG
GTACAAGCGGGCATCCATGCAGGCATCTCACCAGGTCCTGGCAACTCCTACCCCTTGGTGGGACATG
CACCTGGGGCAATTGGTTAGGTGGTGGGCAGGGCCAGGGTGATGGGCTGCAAGGGACTGTCCACAG
GCCCTGCCAGCCTCCACTGAACCCCACTGCTACTTACCCGTCTGACTTTGCCCCCTTAACCCAGC
TCCCTATCCAGCACCCACCTCTTGCAACCAGGTCCTCAGCATGTCTACCCCAAGATGTCTATAAGCCC
AAGGCAGTTGGCAATAGTGGCAACTGGAGGGAGTGGTGCAGGCAAACCTCCCCGGGTGCAGGCACAG
GGGGCAGTGGTGGATATGCAGGCAGTGGGGCAGGGCTTCTACCTGCGACTGCCCAACTGTCAGGAGCT
AGAGCGGCTGGGGCAGCAGCGGCTGGGCTGAGGAAGAAGCCATTACAGCTGCCACATCCCCGGCTGC
GGCAAGGTGTACGGCAAGGCTTCGCATCTGAAAGCCCACTTGCCTGGCACACTGGCGAGAGGCCTTTTCG
TCTGCAACTGGCTTTTCTGTGGCAAGAGTTTACCCGCTCTGACGAGCTGGAGCGCCACGTGCGCACTCA
CACCCGGGAGAAGAAGTTACCTGTCTGCTCTGCTCCAAGCGCTTACCCGAAGCGACCCTTGAGCAAA
CATCAGCGCACCCACGGGAGCCAGGTCCCGGGCCGCCCAAGTGGCCCTAAGGAGCTGGGGGAGGGTC
GCAGCGTCGGGGAAGAAGAAGCCAATCAGCCGCCCGATCTTCCACCTCGCTGCACCCCGAGAAAAGC
CCCTGGAGGCAGCCAGAGCAGCAACCTGCTAGAGATC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



Protein Sequence: >RR200116 representing NM_001037632
 Red=Cloning site Green=Tags(s)


MLTAACSKFGGSSPLRDSTALGKGGTKKPYTDL SAPKTMGDAYPAPF SSTNGLLSPAGSPPAPASGYAND
 YPPFPHSFPGPTGAQDPGLLVKGGHSSSDCLPSVYTSLDMSHPYGSWYKAGIHAGISPGPGNTPTPWDM
 HPGGNWLGQGGQGDGLQGLTSTGPAQPPLNPQLPTYPDFAPLNPAPYPAPHLLQPGPQHVLQPQVYK
 KAVGNSGQLEGSAGKPPRGAGTGGSGGYAGSGAGRSTCDCPNCQELERLGAAAAGLRKKPIHNSCHIPGC
 GKVYGKASHLKAHLRWHTGERPFVCNWLFCGKRFTRSDELERHVRHTTREKKFTCLLCSKRFTRSDHLSK
 HQRTHGEPGPPPPSGPKELGEGRSVGEENQPPRSSTSPAPPEKAPGGSPEQSNLLEI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001037632

ORF Size:	1230 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001037632.1 , NP_001032721.1
RefSeq Size:	1505 bp
RefSeq ORF:	1233 bp
Locus ID:	300260
Cytogenetics:	7q36
MW:	42.8 kDa
Gene Summary:	may act as a transcription factor; mouse homolog is required for bone formation and osteoblast differentiation [RGD, Feb 2006]