

Product datasheet for **RR200126**

Spint2 (NM_001082549) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Spint2 (NM_001082549) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Spint2
Synonyms: HAI-2; Pb
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200126 representing NM_001082549
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGCAGCTATGTGAGCCGCGCGGGCGGGCGCTTCTCGCCCTGTTAGCCTCGCTGCTTCTCTCCG
GGGCCAGGCGGCCAGCAGACCTAGACGTCCACGAGTCTTGTGGAGTTTGAAGGTAGTGGGAAATG
CCGGGCTTCCATCCCAGGTGGTGTACAATGTCCTGATGGGACCTGCCAGCCCTTTGTTATGGAGGC
TGTGAGGGAAATGGTAACAATACTCAAGCAAGGAGGAGTGTCTTGACAAGTGTGCTGGTATCACTGAGA
ACACCATTGATGGTCTGGCCAGGAGCAGTATGGCGCTTCTGCCCTGAGCGTTCCAGAAAAGCAGGATTC
TGCGGACCTCGCTGGTGAGATCTTCAACTATGAAGAATACTGTGTCCAAAGGCGGTCACTGGGCTTGC
CGCGCAGCCTTCCCTCGCTGGTACTATGATGTTGAAAAGAATTCTGTAGCAGCTTATCTATGGCGGGT
GCAGAGGCAACAAGAAGCAGTTACCTCTCCAGGAAGCGTGCATGCAGCACTGCTCTGGGAAGCAGATGTA
TCCTTTCCTGACCCCGGGCACAAAAGCGGTGATCCTGGTGGGGCTGTTTCATCATGGTCTGATCCTCCTC
CTGGGATCCTCCATGGTGTTCCTGATCCGAATGGTTCGCAGAAAGCAGGAGCGTCTCCGCACTGTTT
GGAGCACTGCGGATGACAAGGAGCAGCTGGTGAAGAATACCTGTGTCTTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR200126 representing NM_001082549
 Red=Cloning site Green=Tags(s)

MAQLCEPRRGRALLALLASLLLGAQAASRDLVDHESCGVSKVVGKCRASIPRWYVNTDGTCPFVYGG
 CEGNGNNYPSKEECLDKCAGITENTIDGLARSSDGASALSVPRKQDSADLAGEIFNYEYCVPKAVTGPC
 RAAFPRWYDVEKNSSCFIYGGCRGNKNSYLSQEACMQHCSGKQMYPFLLTPGTKAVILVGLFIMVLILL
 LGSSMVFLIRMVRRKQERALRTVWSTADDKEQLVKNTCVL


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001082549
 ORF Size: 750 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001082549.1 , NP_001076018.1
RefSeq Size:	1401 bp
RefSeq ORF:	753 bp
Locus ID:	292770
Cytogenetics:	1q21
MW:	27.5 kDa