

Product datasheet for RR200149L4

Snx27 (NM_001110151) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Snx27 (NM_001110151) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Snx27
Synonyms:	Mrt1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200149).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 260323
UniProt ID: [Q8K4V4](#)
Cytogenetics: 2q34
Gene Summary: may play a role in neuronal adaptation [RGD, Feb 2006]