

Product datasheet for **RR200151L3V**

Obp3 (NM_001033959) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Obp3 (NM_001033959) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Obp3
Synonyms:	MGC108576
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001033959
ORF Size:	537 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200151).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001033959.2 , NP_001029131.1
RefSeq Size:	873 bp
RefSeq ORF:	540 bp
Locus ID:	259247
Cytogenetics:	5q24
Gene Summary:	odorant binding protein and member of the alpha(2u)-globulin family [RGD, Feb 2006]



[View online »](#)