

Product datasheet for **RR200165L3V**

Fmo4 (NM_144561) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fmo4 (NM_144561) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fmo4
Synonyms:	MGC93155
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_144561
ORF Size:	1680 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200165).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_144561.2 , NP_653147.1
RefSeq Size:	2307 bp
RefSeq ORF:	1683 bp
Locus ID:	246247
UniProt ID:	Q8K4B7
Cytogenetics:	13q22
Gene Summary:	This protein is involved in the oxidative metabolism of a variety of xenobiotics such as drugs and pesticides.[UniProtKB/Swiss-Prot Function]



[View online »](#)