

Product datasheet for **RR200167L3V**

Nexn (NM_139231) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nexn (NM_139231) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nexn
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_139231
ORF Size:	1818 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200167).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_139231.1 , NP_631977.1
RefSeq Size:	1821 bp
RefSeq ORF:	1821 bp
Locus ID:	246172
UniProt ID:	Q9Z2J4
Cytogenetics:	2q45
Gene Summary:	binds actin filaments, b-nexillin isoform also has actin filament cross-linking activity; may play a role in cell- matrix adherens junctions [RGD, Feb 2006]


[View online »](#)