

Product datasheet for **RR200168**

Kcnc2 (NM_139217) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Kcnc2 (NM_139217) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Kcnc2
Synonyms:	KShIIIA; Kv3.2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR200168 representing NM_139217
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGGCAAGATCGAGAACAACGAGAGGGTGATCCTCAATGTCGGAGGCACCAGGCACGAAACCTACCGCA
 GCACTCTCAAGACCCCTTCTGGAACCTCGCCTGGCCCTTCTCGCCTCCTCTGAACCTCAGGGCGACTGCCT
 GACTGCTGCGGGTGACAAGCTGCAGCCGCTGCCCCCTCCGCTGTCTCCACCGCCGCGACCGCCTCCCTTG
 TCCCTGTCCCGAGCGGCTGCTTCGAGGGCGGCGCAGGCAACTGCAGTTCGCACGGTGGCAATGGCAGCG
 ACCACCTGGGGAGGCCGGAATCTTCTTCGATCGCCACCCAGGAGTCTTCGCCTATGTGCTCAACTA
 CTACCGCACGGGCAAGCTGCACTGCCCCCGCAGCTGTGTGGACCGCTTTCGAGGAAGAGCTGGCATT
 TGGGGCATCGATGAGACCGAGCTGGAGCCCTGCTGCTGGATGACCTACAGGCAGCACCGGACGCGGAGG
 AGGCCCTGGATATCTTCGAGACCCGACCTCATCGGAGGCGACCCTGGTATGATGAGGACCTAGGGGG
 CAAGAGACTGGGCAATTGAGGATGCTGCGGGCTGGGAGGACCCGATGGCAAGCTGGCCGCTGGAGGAAG
 CTGCAGCCTCGCATGTGGGCTCTTTGAGGACCCCTATTCATCCAGAGCCGCTAGGTTTATTGCTTTTG
 CTCTCTGTTCTTCAATTTGGTTTCCATCACAACCTTTTGCCTGGAGACACAGAACTTCAATATTGT
 TAAAAACAAGACAGAGCCAGTCAACAGGCACCAGCGCTGTTCTCCAGTATGAAATCGAAACGGATCCT
 GCCTTGACATATGTGGAAGGAGTGTGTGGTGTGGTTACTTTTGAATTTTAGTCCGTATTGTTTTCT
 CGCCCAATAAACTTGAGTTCATCAAAAATCTATTGAACATCATTGACTTTGTGGCCATCCTCCCTTCTA
 CTTAGAGTGAGGACTCAGCGGGCTGTCTCAAAGCGGCTAAAGATGTGCTCGGCTTCTCAGGGTGGTT
 AGGTTTGTGAGGATCCTGAGAATCTTCAAGCTTACCCGCCATTTTCGTAGGTCTGAGAGTCTCGGACACA
 CTCTTCGTGCGAGCACCAATGAATTTTGTGCTGATCATCTTTCTGGCTCTGGGAGTTTGTATTTGCG
 TACGATGATCTACTACGCTGAGCGAGTAGGGGCTCAACCTAATGATCCCTCAGCGAGTGAGCACACACAG
 TTCAAAAACATCCCATTTGGTTTCTGGTGGGCTGTGGTGACCATGACTACCTTAGGCTATGGGGATATG
 ACCCCCAAACATGGTCAGGGATGTTGGTGGGGCCTTGTGTGCTCTGGCTGGAGTCTGACCATAGCTAT
 GCCTGTGCCCGTCATTGTCAACAATTTGGGATGTACTACTCCTTGGCAATGGCGAAGCAGAAACTTCCA
 AGAAAAAGAAAGAAGCACATTCCTGCCCTCTGGCAAGCTCACCTACATTTTGAAGACAGAATTA
 ACATGGCTTGTAACAGTACCCAGAGTGACACATGTCTGGGCAAGAAAACCGGCTTCTGGAACATAACAG
 ATCAGTGTATCAGGTGACGACAGTACAGGAAGTGAAGCCGCTTATCACCTCCGAAAGGCTCCCATC
 AGACGCTCTAGTACCAGAGACAAAAACAGAAGAGGGGAAACATGTTTCTGTTGACGACAGGTGATTACA
 CGTGCGCTTCTGATGGAGGAATCAGGAAAGGATATGAAAAATCCCGAAGCTTAAACAACATAGCGGGCTT
 GGCAGGCAATGCTCTGAGACTCTCTCCAGTAACATCCCCCTACAACCTCTCCGTGTCTCTGAGGGCTCT
 CGGTCTCCCATCCCATCTATCTTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR200168 representing NM_139217
 Red=Cloning site Green=Tags(s)


MGKIENNERVILNVGGTRHETYRSTLTKLTPGTRLALLASSEPGDCLTAAGDKLQPLPPPLSPPRPPPL
 SPVPSGCFEGGAGNCSSHGGNGSDHPGGGREFFDHRPGVFAYVLNYYRTGKLHCPADVCGPLFEEELAF
 WGIDETDVEPCWMTYRQHRDAEEALDIFETPDLIGGDPGDEDLGGKRLGIEDAAGLGGPDGKSGRWRK
 LQPRMWALFEDPYSSRAARFIAFASLFFILVSITTFCLETHEAFNIVKKNKTEPVIINGTSAVLQYEIETDP
 ALTYVEGVCVWFTEFLVRIVFSFNKLEFIKLLNLIIDFVAILPFYLEVGLSGLSSKAADVLGFLRVV
 RFRVILRIFKLRHFVGLRVLGHTLRASNEFLLLIIFLALGVLIFATMIYYAERVAQPNPDSASEHTQ
 FKNIPIGFWWAVVTMTLGYGDMYPQTSWGLVLCALAGVLTAMPVPIVNNFGMYSLAMAKQKLP
 RKRKKHIPPAPLASSPTFCKTELNMACNSTQSDTCLGKENRLLHNRVLSGDDSTGSEPLSPPERLPI
 RRSSTRDKNRGETCFLTTGDYTCASDGGIRKGYEKSRLNLIAGLAGNALRLSPVTSPLYNSPCPLRRS
 RSPIPSIL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_139217

ORF Size: 1914 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_139217.1](#), [NP_631963.1](#)

RefSeq Size:	2441 bp
RefSeq ORF:	1917 bp
Locus ID:	246153
UniProt ID:	P22462
Cytogenetics:	7q22
MW:	70.2 kDa
Gene Summary:	a voltage-gated K ⁺ channel; functions as a delayed rectifier-type K ⁺ channel activated by large membrane depolarizations [RGD, Feb 2006]