

## Product datasheet for **RR200172L4V**

### Zc3h14 (NM\_001033951) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Zc3h14 (NM_001033951) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Zc3h14
Synonyms:	Npuk68
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001033951
ORF Size:	2208 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200172).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001033951.1</a> , <a href="#">NP_001029123.1</a>
RefSeq Size:	3449 bp
RefSeq ORF:	2211 bp
Locus ID:	192359
UniProt ID:	<a href="#">Q7TMD5</a>
Cytogenetics:	6q32
Gene Summary:	recognized by an antibody against a liver cell chromatin fraction [RGD, Feb 2006]



[View online »](#)