

## Product datasheet for **RR200187L3V**

### Reg3a (NM\_001145846) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Reg3a (NM_001145846) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Reg3a
Synonyms:	Pap2; PapII
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001145846
ORF Size:	522 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200187).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001145846.2</a> , <a href="#">NP_001139318.1</a>
RefSeq Size:	805 bp
RefSeq ORF:	525 bp
Locus ID:	171162
UniProt ID:	<a href="#">P35231</a>
Cytogenetics:	4q33
Gene Summary:	lectin-related secretory protein; overexpressed during the acute phase of pancreatitis [RGD, Feb 2006]



[View online »](#)