

## Product datasheet for **RR200189L3V**

### **Lnpep (NM\_001113403) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Lnpep (NM_001113403) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Lnpep
Synonyms:	GP160; IRAP; OTase; P-LAP; Vp165
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001113403
ORF Size:	3075 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200189).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001113403.1</a>
RefSeq Size:	3191 bp
RefSeq ORF:	3078 bp
Locus ID:	171105
UniProt ID:	<a href="#">P97629</a>
Cytogenetics:	1q12
Gene Summary:	a zinc-dependent aminopeptidase that interacts with acyl-coenzyme A (CoA) dehydrogenases and plays a role in retention of GLUT4-containing vesicles to an intracellular compartment [RGD, Feb 2006]



[View online »](#)