

Product datasheet for RR200200L4

Tcf3 (NM_133524) Rat Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tcf3 (NM_133524) Rat Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	Tcf3
Synonyms:	Pan2; Tcfe2a
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200200).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.



[View online »](#)

Locus ID: 171046
UniProt ID: [P21677](#)
Cytogenetics: 7q11
Gene Summary: transcriptional regulator; binds enhancer motifs of pancreatic exocrine genes [RGD, Feb 2006]