

## Product datasheet for **RR200209**

### **Kcnk2 (NM\_172041) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Kcnk2 (NM\_172041) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Kcnk2  
**Synonyms:** rTREK1d; Trek-1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR200209 representing NM\_172041  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGATGAACCCACGGGCGAAAAGGAGCGTCTACCTGGCGGCCCTGACTTGTGGATCCCAAGTCTGCTG  
CTCAGAACTCAAACCGAGGCTCTCGTTCTCCGAAAACCCACCGTCTTGTCCCGGGTGAGAGTGA  
CTCGGCCATTAATGTTATGAAATGGAAGACGGTCTCCACGATTTCTGGTGGTCTCTACCTGATC  
ATCGGAGCCACGGTGTCAAGGCGTTGGAGCAGCCTCAGGAGATTTCTCAGAGGACCACATTGTGATCC  
AGAAACAGAACTTCATAGCCCAGCATGCCTGCGTCAACTCCACCGAGCTGGATGAACTCATCCAGCAAA  
AGTGACGGCCATAAATGCAGGGATTATCCCTTAGGAAACAACCTCCAATCAAGTTAGTCACTGGGACCTC  
GGAAGCTCTTTCTTTGCCGGCACTGTTATCACAACCATAGGATTTGGAAACATCTCCCACGAACTG  
AAGGTGGAAAAATATTCTGTATCATCTATGCCTTGTGGGAATCCCTCTTTGGTTTTCTACTGGCTGG  
GGTTGGGGATCAGCTTGGAAACCATATTTGGAAAAGGAATTGCCAAAGTGGAGGACACATTTAATGAGT  
AATGTTAGTCAGACCAAGATTCGTATCATCTCGACCATCATCTTCATCCTGTTGGCTGTCTCTTCG  
TGGCTCTCCCCGCCGTCATATTCAGCACATAGAAGCTGGAGTGCCTGGACCCATCTACTTTGTGGT  
CATCACTCTGACCACATTGGATTTGGCGATTATGTGGCAGGTGGTCCGACATTGAATATCTGGACTTC  
TACAAGCCGTTGTGTGTTCTGGATCCTCGTTGGGCTGGCCTACTTTGGGCTGTTCTGAGCATGATTG  
GAGACTGGCTACGGTGATATCTAAGAAGACGAAGGAAGAGTGGGAGAGTTCAGAGCGCATGCCGCTGA  
GTGGACAGCCAATGTCACAGCCGAGTTC AAGGAAAACAAGGAGGCGGCTGAGTGTGGAGATCTATGACAAG  
TTCCAGCGTGCCACGTCCGTGAAGCGGAAGCTCTCTGCAGAGCTGGCGGTAACCATAACCAGGAACTGA  
CTCCATGTAGGAGGACCCTGTCGGTGAACCACCTGACCAGCGAGAGGGAAGTCTGCCTCCCTTGTGAA  
GGCTGAGAGCATCTATCTGAACGGTCTGACACCACACTGTCTCGTGAAGACATCGCTGTCATTGAGAAC  
ATGAAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



**Protein Sequence:** >RR200209 representing NM\_172041  
Red=Cloning site Green=Tags(s)


MMNPRAKRSVYLAAPDLLDPKSAQAQNSKPRLSFSKPTVLASRVESDSAINVWKWTVSTIFLVVLYLI  
 IGATVFKALEQPQEISQRTTIVIQKQNFIAQHACVNSTELDELIIQIVTAINAGIPLGNSNQVSHWDL  
 GSSFFFAGTVITIGFGNISPRTEGGKIFCIIYALLGIPLFGLLAGVGDQLGTIFGKGIKVEDTFIKW  
 NVSQTKIRIISTIIIFILFGCVLFVALPAVIFKHIEGWSALDAIYFVVITLTTIGFDYVAGSDIEYLD  
 YKPVVWFVILVGLAYFAAVLSMIGDWRVSKKTKKEEVGEFRAHAAEWTANVTAEFKTRRRLSVEIYDK  
 FQRATSVKRKLSAELAGNHNQELTPCRRTL SVNHLT SEREVL PPLLKAESIYLNGLTPHCAREDI AVIEN  
 MK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:**

NM\_172041

|                               |   |
|-------------------------------|---|
| <b>ORF Size:</b>              | 1266 bp   |
| <b>OTI Disclaimer:</b>        | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>  |
| <b>OTI Annotation:</b>        | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <a href="#">NM_172041.2</a> , <a href="#">NP_742038.2</a>   |
| <b>RefSeq Size:</b>           | 1371 bp   |
| <b>RefSeq ORF:</b>            | 1269 bp   |
| <b>Locus ID:</b>              | 170899  |
| <b>UniProt ID:</b>            | <a href="#">Q920B6</a>  |
| <b>Cytogenetics:</b>          | 13q26   |
| <b>MW:</b>                    | 46.8 kDa  |
| <b>Gene Summary:</b>          | a large conductance K+ channel [RGD, Feb 2006]  |