

Product datasheet for **RR200217L3V**

Tpm3 (NM_057208) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Tpm3 (NM_057208) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Tpm3
Synonyms:	TM30nm; Tpm5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_057208
ORF Size:	741 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200217).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_057208.2 , NP_476556.2
RefSeq Size:	1078 bp
RefSeq ORF:	744 bp
Locus ID:	117557
UniProt ID:	Q63610
Cytogenetics:	2q34
Gene Summary:	rod-like protein associated with actin filament in nonmuscle cells; may mediate the stability of actin filaments; may play a role in the regulation of cell transformation, cytokinesis, motility, and morphogenesis [RGD, Feb 2006]



[View online »](#)