

Product datasheet for **RR200222L3V**

Slc15a1 (NM_001079838) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc15a1 (NM_001079838) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc15a1
Synonyms:	Pept1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001079838
ORF Size:	450 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200222).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001079838.1 , NP_001073307.1
RefSeq Size:	626 bp
RefSeq ORF:	453 bp
Locus ID:	117261
UniProt ID:	P51574
Cytogenetics:	15q24
Gene Summary:	acts as a proton dependent peptide transporter; may transport oral beta-lactam antibiotics across intestinal brush borders [RGD, Feb 2006]



[View online »](#)