

## Product datasheet for **RR200229L4V**

### Calcr (NM\_001034015) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Calcr (NM_001034015) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Calcr
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001034015
ORF Size:	1488 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200229).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001034015.1</a> , <a href="#">NP_001029187.1</a>
RefSeq Size:	3433 bp
RefSeq ORF:	1491 bp
Locus ID:	116506
UniProt ID:	<a href="#">P32214</a>
Cytogenetics:	4q13
Gene Summary:	receptor for calcitonin; acts via adenylate cyclase signaling [RGD, Feb 2006]



[View online »](#)