

## Product datasheet for **RR200236L3V**

### Ugt1a6 (NM\_057105) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | Lentiviral Particles                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Product Name:             | Ugt1a6 (NM_057105) Rat Tagged ORF Clone Lentiviral Particle                                                                                                                                                                                                                                                                                                                                                                                        |
| Symbol:                   | Ugt1a6                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Synonyms:                 | Udpgt; UDPGT 1-6; Ugt1; Ugt1a7                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Mammalian Cell Selection: | Puromycin                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Vector:                   | pLenti-C-Myc-DDK-P2A-Puro (PS100092)                                                                                                                                                                                                                                                                                                                                                                                                               |
| Tag:                      | Myc-DDK                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ACCN:                     | NM_057105                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| ORF Size:                 | 1590 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RR200236).                                                                                                                                                                                                                                                                                                                                                                                     |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.                                                                                                                                                                                                                                                                                                               |
| RefSeq:                   | <a href="#">NM_057105.3</a> , <a href="#">NP_476446.2</a>                                                                                                                                                                                                                                                                                                                                                                                          |
| RefSeq Size:              | 2313 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| RefSeq ORF:               | 1593 bp                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Locus ID:                 | 113992                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Cytogenetics:             | 9q35                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Gene Summary:             | catalyzes the glucuronidation of 4-nitrophenol [RGD, Feb 2006]                                                                                                                                                                                                                                                                                                                                                                                     |


[View online »](#)