

Product datasheet for **RR200237**

Sqstm1 (NM_175843) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Sqstm1 (NM_175843) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sqstm1
Synonyms:	Osi; ZIP; ZIP3
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>NM_175843 ORF sequence, RR200237 may differ due to SNPs. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCTTCGCTCACGGTGAAGGCCTATCTACTGGGCAAGGAGGAGGCGGCCCGGAGATCCGCCGCTTC
AGCTTCTGCTTCAGCCCGGAGCCGGAGGCGGAAGCCGGCTGGCCCGGGCCCTCGCAGAGGCTGCTG
AGCCGGGTGGCTGTGCTGTTCCCGCGCTGCGGCTGGAGGCTTTCAGGGCAGCTACCGCATGAGGAT
GGGGACTTGGTCGCTTCTCCAGTGATGAGGAAGTACAATGGCCATGTCTATGTGAAAGATGACATC
TTCCGCATCTACATTAAGAGAAGAAGGAGTGCCGGCGGGAACATCGCCCCCATGTGCTCAGGAGGCA
CGAAGCATGGTGCACCCAACGTGATTTGTGATGGTTGCAATGGGCTGTGGTGGGAAGTGCCTATAAG
TGCAGTGTGTGCCCGACTACGACCTGTGCAGCGTCTGCGAGGGGAAGGGCTGCACAGGGAGCACAGC
AAGCTCATCTTTCCCAACCCCTTTGGCCACCTCTGTATAGCTTCTCTCATAGCCGCTGGCTTGGGAAG
CTGAAACATGGGCACTTTGGCTGGCCTGGCTGGGAGATGGGCCACCAGGGAAGTGGAGCCACGTCCT
CCTCGCGCAGGGGATGGTCGCCCTTGGCCACAGCTGAGTCGGCTTCTGCTCCATCAGAGGATCCCAAT
GTCAATTTCTGAAGAAATGTGGGGAGAGCGTGGCAGCTGCCCTCAGCCCTTAGGCATCGAGGTTGAC
ATTGATGTGGAACATGGAGGGAAGAGAAGCCGCTGACACCCACTCTGCAGAAAGTTCAGCACAGGC
ACAGAAGATAAGAGTGGTACTCAGCCAAGCAGCTGCTCTTCGGAAGTCAAGCAACCTGACGGGGCCGGG
GAGGGCCCTGCTCAGTCTCTGACAGAGCAGATGAAGAAGATAGCCTTGGAGTCGGTGGGACAGCCAGAG
GAACTGATGGAGTCGGATAACTGCTCAGGAGGGGATGACGACTGGACGATTTGTCTTCTAAAGAAGTG
GACCCATCCACAGGTGAACTCCAGTCTCTACAGATGCCAGAATCGGAAGGGCCAAGCTCTCTAGACCCC
TCACAGGAAGGCCCCACAGGGCTGAAGGAAGCTGCCCTGTACCCACATCTCCACCAGAGGCTGATCCC
CGGCTGATTGAGTCACTCTCTCAGATGCTGTCATGGGTTTCTCGGATGAAGGCGGCTGGCTCACCAGG
CTCCTACAGACCAAGAATTATGACATCGGGGCTGCTCTGGACACGATCCAGTATTCAAAGCACCCCTCA
CCATTG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```




Protein Sequence: >Peptide sequence encoded by RR200237
 Blue=ORF Red=Cloning site Green=Tag(s)


MASLTVKAYLLGKEEAAREIRRFSCFSPEPEAEAAAGPGPCERLLSRVAVLFPALRPGGFQAHYRDED
 GDLVAFSSDEELTMAMSYVKDDIFRIYIKEKKECRREHRPPCAQEARSVMHPNVICDGCNGPVVGYTRYK
 CSVCPDYDLCSVCEGKGLHREHSLIFPNPFGLHSDSFHSRWLRKLRKHGHFGWPGWEMGPPGNWSPRP
 PRAGDGRPCPTAESASAPSEDPNVNFKNVGESVAAAALSPLEGIEVDIDVEHGGKRSRLTPTSAESSSTG
 TEDKSGTQPSSCSSEVSKPDGAGEGPAQSLTEQMKKIALESVGQPEELMESDNCSSGGDDWTHLSSKEV
 DPSTGELQSLQMPSESGPSSLDPSQEGPTGLKEAALYPHLPPEADPRLIESLSQMLSMGFSDGGWLTR
 LLQTKNYDIGAALDTIQYSKHPPPL
 TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_175843

ORF Size:	1317 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_175843.5
RefSeq Size:	1991 bp
RefSeq ORF:	1320 bp
Locus ID:	113894
UniProt ID:	O08623
Cytogenetics:	10q22
MW:	47.7 kDa
Gene Summary:	interacts with PKC-zeta and GABA(C) receptors to form ternary postsynaptic complex in the retina [RGD, Feb 2006]