

Product datasheet for **RR200247**

Slc44a1 (NM_053492) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc44a1 (NM_053492) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc44a1
Synonyms:	Cdw92; Ctl1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR200247 representing NM_053492
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGGATCGCC

ATGGGCTGCTGCAGCTCCGCCTCCGCCGCGCAGAGCTCCAACGAGAATGGAAGCCGTTGGAGGACCGTA
GCTGCACAGACATACCATGGCTGCTGCTTTCGTTCTGTTCTGCATTGGGATGGGATTTATTTGGGCTT
TTCAGTAGCAACAGGTGCAGCAGCAAGACTAGTGTGAGGATACGACAGCTATGGAAATATCTGTGGGCAG
AGAAATGCAAACTGGAAGCAATAGCGAACAGTGGCCTGGACCACACCCACCGGAAGTATGTGTTCTTTT
TGGACCATGCAATCTGGACTTGATAAACAGGAAGATCAAGTCTATGGCACTGTGTAGCGGCATGTCC
ACGACAGGAGTTGAAAACCTGAGTGACGTCCAGAAGTTTGCCGAGATCAATGGTTCAGCACTGTGTAGC
TACAACATAAAGCCTCCGAGTACACTGACTGCAAAATCCTCTGCCTTTGCCCAAACCTACCACTAC
CAGCCAGTGCACCTATCCATTCTCCATCGCTGTGCTCCTGTGAACATTTCTGCTATGCCAAGTTTGC
AGAGGCCCTGATCACCTTTGTCAGTGACAATAGTGTCTTACACAGGCTGATTAGTGGAGTAAAGACCAGC
AAAGAAATATATTGGGACTTTGCTTGTATCACTAGTTCTGTCCATGATTTTGTATGGTATTATCCGGT
ATATCTCAAGAGTGTAGTGTGGATCCTAACGATTCTGGTTCATCCTGGGTTCACTAGGTGGCACCCGGT
ACTATGGTGGCTATATGCAAAGCAGAGAAGTCCCCAAAGAAACAGTTATTCCTGAACAGCTTCAGATC
GCGGAAGACAATCTGCGGGCGCTCCTCATCTATGCCATTTCTGCTACGGTGTTCAGTGTCTCTCTTCC
TGATCATGCTGGTTATGCGTAAACGCGTCGCTCTACCATCGCCTTGTCCACGTGGCTGGCAAGGTCTT
CATACACCTGCCTCTGCTAGTCTCCAGCCCTTCTGGACTTTTTTGTCTCTGCTTATTTTGGGCATAC
TGGATAATGACACTTCTCTTCTTGGCACTACTGGCAGTGTTCAGAATGAACAAGGTTTCGTGGAGT
ACAAAATCTCTGGGCCTTTGCAGTACATGTGGTGGTACCACGTGGTCCGCTGATTTGGATCAGTGAGTT
CATTCTCGCCTGTCAGCAGATGACTGTGGCAGGAGCTGTGGTAACTACTACTTTACAAGGGATAAACGG
AATTTGCCTTTACACCTATTCTGGCATCAGTGAACCGCCTTATTCGTTACCACCTCGGTAAGTGTGGCA
AAGGATCCTTCATTATCACCTCGTCAAGATCCCAGGATGATCCTTATGTATATCCACAGTCAAGTCAA
AGGGAAGGAAAACGCTTGTGCACGATGCATGCTGAAATCCTGCATTTGCTGCCTTTGGTGTCTTAAAAAG
TGCTGAGCTACTTAAACCAGAATGCATACACAGCCACAGCCATCAATAGCACCAACTTCTGCACCTCAG
CAAAGGATGCCTTCGTCATTCTAGTGGAGAATGCTCTGCGAGTGGCTGCCATCAACACAGTGGGGGACTT
CATGTTATTCCTAGGCAAGGTGCTGATAGTCTGCAGCACGGCCTGGCTGGGATTATGCTGCTCAACTAC
CAGCAAGATTACACAGTGTGGGTGCTGCCTCTCATCATCGTCTGCCTTTTGTCTTCTAGTCGCTCATT
GCTTCTGTCCATTTACGAAATGGTAGTGGATGTGTTGTTTTGTGTTTTGCCATTGATACAAAATACAA
TGATGGGAGCCCAGGCCGAGAATTCTATATGGACAAAGTCTAATGGAGTTTGTGGAAAACAGTAGGAAA
GCAATGAAAGAAGCTGTTAAGGGAGGCCCGCGGATGCCCGAGAGCTAAAGCCGATGCTCAGGAAAAGG

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR200247 representing NM_053492
 Red=Cloning site Green=Tags(s)


MGCCSSASAAQSSKREWKPLEDRSCTDIPWLLLFVLFICIGMGFICGFSVATGAAARLVSGYDSYGNICGQ
 RNAKLEAIANGLDHTHRKYVFFLDPCNLDLNRKIKSMALCVAACPRQELKTLSDVQKFAEINGSALCS
 YNIKPSEYTLTAKSSAFCPKLPVPASAPIPFHRCAPVNI SCYAKFAEALITFVSDNSVLHRLISGVMTS
 KEIILGLCLLSLVL SMILMVIIRYISRVLVWILTILVILGSLGGTGVLWVLYAKQRRSPKETVIPEQLQI
 AEDNLRALLIYAI SATVFTVILFLIMLVMRKRVALTIALFHVAGKVF IHLPLLVFQPFWTFALVLFWAY
 WIMTLLFLGTTGSAVQNEQGFVEYKISGPLQYMWWYHVVL IWI SEF ILACQMQMTVAGAVVYYFTRDKR
 NLPFTPILASVNRLIRYHLGTVAKGSFIITLVKIPRMILMYIHSQKKGKENACARCMLKSCICCLWCLEK
 CLSYLNQAYTATAINSTNFCTSAKDAFVILVENALRVA AINTVGDFMLFLGKVLIVCSTGLAGIMLLNY
 QQDYTVVWVLP LIIVCLFAFLVAHCFLSIYEMVVDVLF LCF AIDTKYNDGSPGREFYMDKVLMEFVENS RK
 AMKEAGKGAADARELKPMLRKR

SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_053492

ORF Size: 1959 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_053492.3](#), [NP_445944.2](#)

RefSeq Size:	2851 bp
RefSeq ORF:	1962 bp
Locus ID:	85254
UniProt ID:	<u>Q8VII6</u>
Cytogenetics:	5q24
MW:	73.1 kDa
Gene Summary:	may be a choline transporter for a specific subset of cholinergic neurons [RGD, Feb 2006]