

## Product datasheet for **RR200247L3V**

### Slc44a1 (NM\_053492) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc44a1 (NM_053492) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc44a1
Synonyms:	Cdw92; Ctl1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053492
ORF Size:	1959 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200247).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053492.3</a> , <a href="#">NP_445944.2</a>
RefSeq Size:	2851 bp
RefSeq ORF:	1962 bp
Locus ID:	85254
UniProt ID:	<a href="#">Q8VII6</a>
Cytogenetics:	5q24
Gene Summary:	may be a choline transporter for a specific subset of cholinergic neurons [RGD, Feb 2006]



[View online »](#)