

Product datasheet for **RR200248**

Slc44a1 (NM_001033852) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc44a1 (NM_001033852) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc44a1
Synonyms:	Cdw92; Ctl1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>RR200248 representing NM_001033852
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGGATCGCC

ATGGGCTGCTGCAGCTCCGCCTCCGCCGCGCAGAGCTCCAACGAGAATGGAAGCCGTTGGAGGACCGTA
 GCTGCACAGACATACCATGGCTGCTGCTTTCGTTCTGTTCTGCATTGGGATGGGATTTATTTGTGGCTT
 TTCAGTAGCAACAGGTGCAGCAGCAAGACTAGTGTGAGGATACGACAGCTATGGAATATCTGTGGGCGAG
 AGAAATGCAAACTGGAAGCAATAGCGAACAGTGGCCTGGACCACCCCACCGGAAGTATGTGTTCTTTT
 TGGACCATGCAATCTGGACTTGATAAACAGGAAGATCAAGTCTATGGCACTGTGTAGCGGCATGTCC
 ACGACAGGAGTTGAAAACCTGAGTGACGTCCAGAAGTTTGCCGAGATCAATGGTTCAGCACTGTGTAGC
 TACAACATAAAGCCTCCGAGTACACTCTGACTGCAAAATCCTCTGCCTTTGCCCAAACCTACCACTAC
 CAGCCAGTGCACCTATCCATTCTCCATCGCTGTGCTCTGTGAACATTTCTGCTATGCCAAGTTTGC
 AGAGGCCCTGATCACCTTTGTCAGTGACAATAGTGTCTTACACAGGCTGATTAGTGGAGTAAAGACCAGC
 AAAGAAATATATTGGGACTTTGCTTGTATCACTAGTTCTGTCCATGATTTTGTATGGTATTATCCGGT
 ATATCTCAAGAGTGTAGTGTGGATCCTAACGATTCTGGTTCATCTGGGTTCACTAGGTGGCACCCTGT
 ACTATGGTGGCTATATGCAAAGCAGAGAAGTCCCCAAAGAAACAGTTATTCCTGAACAGCTTCAGATC
 GCGGAAGACAATCTGCGGGCGCTCCTCATCTATGCCATTTCTGCTACGGTGTTCAGTGTCTCTCTTCC
 TGATCATGCTGGTTATGCGTAAACGCGTCTGCTCACCATCGCCTTGTCCACGTTGGCTGGCAAGTCTT
 CATACACCTGCCTCTGCTAGTCTCCAGCCCTTCTGGACTTTTTTGTCTCTGCTTATTTTGGGCATAC
 TGGATAATGACACTTCTCTTCTTGGCACTACTGGCAGTGTTCAGAATGAACAAGTTTCGTGGAGT
 ACAAATCTCTGGGCCTTTGCAGTACATGTGGTGGTACCACGTGGTCCGCTGATTTGGATCAGTGAGTT
 CATTCTCGCCTGTCAGCAGATGACTGTGGCAGGAGCTGTGGTAACTACTACTTTACAAGGATAAACGG
 AATTTGCCTTTACACCTATTCTGGCATCAGTGAACCGCCTTATTCGTTACCACCTCGGTACTGTGGCCA
 AAGGATCCTTCATTATCACCTCGTCAAGATCCCAGGATGATCCTTATGTATATCCACAGTCAAGTCAA
 AGGGAAGGAAAACGCTTGTGCACGATGCATGCTGAAATCCTGCATTTGCTGCCTTTGGTGTCTTAAAAAG
 TGCTGAGCTACTTAAACCAGAATGCATACACAGCCACAGCCATCAATAGCACCAACTTCTGCACCTCAG
 CAAAGGATGCCTTCGTCATTCTAGTGGAGAATGCTCTGCGAGTGGCTGCCATCAACACAGTGGGGGACTT
 CATGTTATTCCTAGGCAAGGTGCTGATAGTCTGCAGCACGGCCTGGCTGGGATTATGCTGCTCAACTAC
 CAGCAAGATTACACAGTGTGGGTGCTGCCTCTCATCATGCTCTGCCTTTTGTCTTCTAGTCTGCTCATT
 GCTTCTGTCCATTTACGAAATGGTAGTGGATGTGTTGTTTTGTGTTTTGCCATTGATACAAAATACAA
 TGATGGGAGCCCAGGCCGAGAATTCTATATGGACAAAGTCTAATGGAGTTTGTGGAAAACAGTAGGAAA
 GCAATGAAAGAAGCTGTTAAGGGAGGCCCGCGGATGCCGAGAGCTAAAGCCGATGGCTTCGGGAGCAA
 GTTCTGCT

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR200248 representing NM_001033852
 Red=Cloning site Green=Tags(s)


MGCCSSASAAQSSKREWKPLEDRSCTDIPWLLLFVLFICIGMGFICGFSVATGAAARLVSGYDSYGNICGQ
 RNAKLEAIANGLDHTHRKYVFFLDPCNLDL INRRIKSMALCVAACPRQELKTLSDVQKFAEINGSALCS
 YNIKPSEYTLTAKSSAFCKLPVPASAPIPFHRCAPVNI SCYAKFAEALITFVSDNSVLHRLISGVMTS
 KEIILGLCLLSLVL SMILMVIIRYISRVLVWILTILVILGSLGGTGVLWVLYAKQRRSPKETVIPEQLQI
 AEDNLRALLIYAI SATVFTVILFLIMLVMRKRVALTIALFHVAGKVF IHLPLLVFQPFWTFALVLFWAY
 WIMTLLFLGTTGSAVQNEQGFVEYKISGPLQYMWYHVVGLIWI SEFILACQMQTVAGAVVTYFTRDKR
 NLPFTPI LASVNRLIRYHLGTVAKGSFIITLVKIPRMILMYIHSQKKGKENACARCMLKSCICCLWCLEK
 CLSYLNQAYTATAINSTNFCTSAKDAFVILVENALRVA AINTVGFMLFLGKVLIVCSTGLAGIMLLNY
 QQDYTVVWLPLIIVCLFAFLVAHCFLSIYEMVVDVFLCFAIDTKYNDGSPGREFYMDKVLMEFVENSRK
 AMKEAGKGAADARELKPMASGASSA

SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001033852

ORF Size: 1968 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001033852.1](#), [NP_001029024.1](#)

RefSeq Size:	4507 bp
RefSeq ORF:	1971 bp
Locus ID:	85254
UniProt ID:	<u>Q8VII6</u>
Cytogenetics:	5q24
MW:	73.1 kDa
Gene Summary:	may be a choline transporter for a specific subset of cholinergic neurons [RGD, Feb 2006]