

## Product datasheet for **RR200252L4V**

### Trpc4 (NM\_080396) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Trpc4 (NM_080396) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Trpc4
Synonyms:	Trp4; Trrp4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_080396
ORF Size:	2931 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200252).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_080396.1</a> , <a href="#">NP_536321.1</a>
RefSeq Size:	3127 bp
RefSeq ORF:	2934 bp
Locus ID:	84494
Cytogenetics:	2q26
Gene Summary:	forms a receptor regulated cation channel; mediates calcium transport in response to histamine H1 receptor activation [RGD, Feb 2006]



[View online »](#)