

## Product datasheet for **RR200274L2V**

### Cdk16 (NM\_031077) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdk16 (NM_031077) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cdk16
Synonyms:	Pctk1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_031077
ORF Size:	1383 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200274).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031077.1</a> , <a href="#">NP_112339.1</a>
RefSeq Size:	3111 bp
RefSeq ORF:	1386 bp
Locus ID:	81741
UniProt ID:	<a href="#">Q63686</a>
Cytogenetics:	Xq11
Gene Summary:	cyclin-dependent kinase related protein that may function in differentiation [RGD, Feb 2006]



[View online »](#)