

Product datasheet for **RR200279L1V**

Ceacam1 (NM_001033860) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ceacam1 (NM_001033860) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ceacam1
Synonyms:	Bgp; BGPR; Ccam1; CD66a
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_001033860
ORF Size:	1557 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200279).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001033860.1 , NP_001029032.1
RefSeq Size:	3573 bp
RefSeq ORF:	1560 bp
Locus ID:	81613
UniProt ID:	P16573
Cytogenetics:	1q21
Gene Summary:	glycoprotein component of bile; plays a role in cell adhesion [RGD, Feb 2006]



[View online »](#)