

Product datasheet for **RR200297L3V**

Hnrnpd (NM_001082541) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Hnrnpd (NM_001082541) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Hnrnpd
Synonyms:	Auf1; Hnrpd
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001082541
ORF Size:	855 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200297).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001082541.1 , NP_001076010.1
RefSeq Size:	1514 bp
RefSeq ORF:	858 bp
Locus ID:	79256
UniProt ID:	Q9J154
Cytogenetics:	14p22
Gene Summary:	mouse homolog is an RNA binding protein that binds A+U-rich elements [RGD, Feb 2006]



[View online »](#)