

Product datasheet for **RR200312**

Bcl2l11 (NM_171989) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Bcl2l11 (NM_171989) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Bcl2l11
Synonyms: Bim; BimL; Bod
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200312 representing NM_171989
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

**ATGGCCAAGCAACCTTCTGATGTAATTCTGAGTGTGACAGAGAAGGTGGACAATTGCAGCCTGCTGAGA
GGCCTCCCAGCTCAGGCCTGGGGCCCTACCTCCCTACAGACAGAATCGCAAGCTTCATAAGGCAGTC
TCAGGAGGACCTGAAGATTTGCGCCAGAGATACGGATCGCACAGGAGCTGCGGCGGATCGGAGACGAG
TTCAATGAGACTTACACGAGGAGGGCGTTTGCAAACGATTACCGAGAGCGGAAGACCACCCGAAATGG
TTATCTTACAACCTGTTACGATTCATCTCCGTCTGGTCTGGAGAAGGCAC**

**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA**

Protein Sequence: >RR200312 representing NM_171989
Red=Cloning site Green=Tags(s)

MAKQPSDVNSECDREGGQLQPAERPPQLRPGAPTSLQTESQASIRQSQEDPEDLRPEIRIAQELRRIGDE
FNETYTRRAFANDYREAEDHPQMVLQLLRFIFRLVWRRH


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: Sgfl-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_171989

ORF Size: 330 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_171989.1</u> , <u>NP_741986.1</u>
RefSeq Size:	333 bp
RefSeq ORF:	333 bp
Locus ID:	64547
UniProt ID:	<u>O88498</u>
Cytogenetics:	3q36
MW:	12.9 kDa
Gene Summary:	proapoptotic Bcl-2 related protein that interacts with anti-apoptotic Bcl-2 family members to play a role in regulation of apoptosis [RGD, Feb 2006]