

## Product datasheet for **RR200313L3V**

### **Bcl2l11 (NM\_022612) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Bcl2l11 (NM_022612) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Bcl2l11
Synonyms:	Bim; BimL; Bod
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022612
ORF Size:	420 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR200313).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022612.1</a> , <a href="#">NP_072134.1</a>
RefSeq Size:	423 bp
RefSeq ORF:	423 bp
Locus ID:	64547
UniProt ID:	<a href="#">O88498</a>
Cytogenetics:	3q36
Gene Summary:	proapoptotic Bcl-2 related protein that interacts with anti-apoptotic Bcl-2 family members to play a role in regulation of apoptosis [RGD, Feb 2006]



[View online »](#)