

Product datasheet for **RR200331**

Ensa (NM_001033974) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ensa (NM_001033974) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ensa
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR200331 representing NM_001033974
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCAGAAACAAGAAGAAGAAAACCTGCGGAGGAAACCGGCGAGGAGAAGCAGGATACACAGGAGA
AAGAAGGTATTCTTCTGAGAAAGCTGAGGAAGCAAAGCTAAAGGCCAAATATCCAAGCCTAGGACAAAA
GCCTGGAGGCTCCGACTTCCTCATGAAGAGACTCCAGAAAGGGCAAAGTACTTTGACTCAGGAGACTAC
AACATGGCCAAAGCCAAGATGAAGAACAAGCAGCTACCGAGTGCAGGAGCAGACAAGAACCTGGTGACCG
GTGACCACATCCCCACCCACAGGACCTGCCCCAGAGAAAGTCTCGCTCGTCACCAGCAAGCTTGCGGG
G

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR200331 representing NM_001033974
Red=Cloning site Green=Tags(s)

MSQKQEEENPAEETGEEKQDTQEKEGILPEKAEAKLKAKYPSLGQKPGGSDFLMKRLQKQKYFDSDGY
NMAKAKMKNKQLPSAGADKNLVTGDHIPTQDLPQRKSSLVTSKLAG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001033974

ORF Size: 351 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001033974.2</u> , <u>NP_001029146.1</u>
RefSeq Size:	2296 bp
RefSeq ORF:	354 bp
Locus ID:	60334
UniProt ID:	<u>P60841</u>
Cytogenetics:	2q34
MW:	12.9 kDa
Gene Summary:	may mediate insulin secretion through interaction with the pancreatic beta-cell ATP-sensitive potassium (K(ATP)) channel [RGD, Feb 2006]